

8.0 PARTS LIST

Description	Part Number
Filter 0.2 Micron 2.5 X 10	20-3021
Flow Control, 2.0GPMX.5MPT	42-0043
Fuse, 10 amp, FNM	RFNQ-R-10
Fuse, 2A	63760133
Hose Assy, Bicarb, Inlet, Disinfection	0110112-0103C
Logic Relay Controller	64-0026
Multifunction Digital Timer	64-0025
Peracidin Disinfectant, 2 Quarts	95-0006
Peracidin Disinfectant, 4 Quarts	95-0007
Pumps, .75 HP,1 Phase, 1" Inlet, 1" Outlet	80-0264
Q-Con, CPC, Female Coup x .38MPT, Polysulfone	16-0077
Q-Con, CPC, Female Coup, NSHD17008, With Shut Off, .5 HB	16-0076
Q-Con, CPC, Male Coup Insert, NSHD24008, With Shut Off, .5 MPT	16-0075
Q-Con, CPC, Male Insert x .38MPT, No Shut Off, Polysulfone	16-0078
Relay, 2-POLE, 25 AMP	64-0046
Test Strips Peracetic Acid Test	97HP20401
Test Strips Renal Check, Residual Peroxide	97PX20501
Test Strips Water Check 2 Low Level Chlorine/Chloromine	97CM20201
Test Strips Water Check Residual Chlorine/Chloromine	97RC22101
Test Strips, Ozone (in Water), 0.05 to >0.5 ppm	97K100-0111
Tubing .25OD, Blue	08760287
Tubing, .5X.063, Poly, Blue	08760280
Tubing, 3/4, PT2401BK, Black	08-0007
VAL, Ball, .5, True Union, PVC80, Panel Mount	041530953
VAL, Ball, .75, True Union, PVC80, Panel Mount	041530844
VAL, Check, .5FPT, 1/3#, White, PP, SMC (Return Loop)	55-0003
VAL, Check, .5FPT, PP, 3# Cracking (RO Inlet Water)	55-0010
VAL, SOL, .5FPT, NC, 120V, DIN	59-0002
Valve Ball 3/8" M X F (Sample Port & Drain)	041001
Valve, Ball, EL, .25MPTXFPT, PVC80 (Sample Port)	041004
Valve, Needle, .5, FPT, PVC80	04100002
KIT,BICARB SPRAY NOZZLE ASSY,55 GAL,MIX/DIST	911-9118
KIT,BICARB SPRAY NOZZLE ASSY,100-185 GAL.,MIX/DIST	911-9119
KIT,BICARB,MANIFOLD-PANEL,FEED	911-9120
KIT,BICARB TANK INTERNAL,MIX,55 GALLON	911-9121
KIT,BICARB TANK INTERNAL,MIX,100 GALLON	911-9122