


REV	DATE	DESCRIPTION	INIT
A	07JAN09	ORIGINAL ISSUE	JJJ



047-0001	22	ADPT, FIGHT X 0.5 FPT, PP
14520417	21	ADPT, 0.75 MCHT X 0.5 HB, PP
046520423	20	ADPT, MCHT X 0.5 MPT, PP
14570110	19	ADPT, GH, SN, INS, 0.5 HB, PP
14570111	18	ADPT, GH WASHER, RUBBER
14-0014	17	ADPT, GH SWIVEL NUT, PP
15-0004	16	FERRULE HOSE CLAMP, 0.781" ID
15650050	15	CLAMP, HOSE, 0.5, SST
19530901	14	SILEX AIR VENT
19530908	13	SILEX KNOB
92530122	12	STUD, 5/16 x 3.5 THD, SS
14520403	11	ADPT, EL, .25MPT X .5HB, PP
041002	10	VAL, BALL, .25, MPT X FPT, PVC80
19300001-00	9	O-RING, SILEX, MOLDED UNIT
14520409	8	ADPT., EL, 0.75 MPT X HB, PP
14520406	7	ADPT., EL, .05, MPT X HB, PP
75679105	6	RESISTIVITY LIGHT 10K
94560411	5	CASTER, TWIN WHEEL, NO BRAKE
94560410	4	CASTER, TWIN WHEEL WITH BRAKE
41-0002	3	FLOWMETER, 0.5-5, INLINE
80-0014	2	HWAKI PUMP WMD-30RZT, 115V
01119-0010	1	COLUMN,SILEX,MOLDED ASSY, SINGLE OR SERIES
AWI P/N	ITEM #	DESCRIPTION

REPLACEMENT COMPONENTS

THIS IS A LISTING OF THE MOST COMMON PARTS REPLACED FOR THE 00M10800, 00M20800 & 00M30800 SYSTEMS