

## 7.0 PARTS LIST

Description	Part Number
FILTER CARTRIDGE SUBMICRON, 0.2 MICRON, 2.5" X 10", DOE	20-3021
FLOATSWITCH,ASSY	0167-0010
FUUSE 2A	63760133
FUUSE 15A	63760186
GAUGE,0-30,.25,CBM	43-0006
LIGHT AMBER PILOT 120V	65511227
HORN PIEZO 120 VAC, 80-95 DB	65760204
LIGHT RED PILOT 120V	65511223
LINE CORD 120V 12/3 HOSPITAL GRADE STRIPPED ENDS	66932113
MIX & DISTRIBUTION PUMP	80-0102
PERACIDIN DISINFECTANT, 2 QUARTS	95-0006
PERACIDIN DISINFECTANT, 4 QUARTS	95-0007
RELAY,120V,15AMP,3PDT,11TERM-LED&TEST BUTTON	64-0017
SWITCH 3-POSITION MAINTAIN	65511222
SWITCH 3-POSITION RETURN RIGHT	65511229
SWITCH PUSH BUTTON	65511220
TEST STRIPS PERACETIC ACID TEST	97HP20401
TEST STRIPS RENAL CHECK, RESIDUAL PEROXIDE	97PX20501
TEST STRIPS WATER CHECK RESIDUAL CHLORINE/CHLORAMINE	97RC22101
TEST STRIPS WATER CHECK 2 LOW LEVEL CHLORINE/CHLORAMINE	97CM20201
TIMER RAIL MOUNT 11 PIN	64760235
TRANSFORMER 120/208/240-24V 40VA	62760177
VAL,CHECK,.5FPT,WHITE,PP	55-0003
VAL,SOL,.5FPT,NC,120V,DIN	59-0002
VALVE BALL 1.25" TRUE UNION PVC80	041732001
VALVE BALL 1/2" TRUE UNION PVC80	041531812
VALVE BALL 3/4" TRUE UNION PVC80	041530841
VALVE BALL 3/8" M X F (SAMPLE PORT/DRAIN)	041001
VALVE,NEEDLE,.5,FPT,PVC80 (THROTTLE VALVE)	04100002
WIRE HARNESS,DIN,STD ISO,6FT,STRIPPED ENDS	66932109
WRENCH,FILTER,10" OR 20" WHITE	21675184