













# SECTION 10 - MRO1 SPARE PARTS LIST



NO.	DESCRIPTION	PART #
1	FLOW METER	41530606
2	CONTROLLER ASSEMBLY	R69-0019
3	TEST PORT	041004
4	FLUSH VALVE	04-5001
5	PAA CONNECTOR	16-0042
6	3/8 INCH BULKHEAD	10-L277
7	1/2 INCH BULKHEAD	10-L278
8	PAA BOTTLE ASSEMBLY	0185-0020
9	PUMP PRESSURE GAUGE	430001
10	PREFILTER GAUGES	430002
11	FILTER HOUSING	21-0032
12	FILTER O-RING KIT	0013-0001
13	RO PUMP ASSEMBLY	R080-0002
14	PT401 PUMP (OPTIONAL)	R84-0002
15	MEMBRANE ONLY	R22-4026
16	NEPHROS DSU 510K FILTER	20-0009
17	PRESSURE SWITCH	R65-0014
18	MOTOR STARTER	R61-0008
19	PRESSURE REGULATOR	R44-0066
20	CONDUCTIVITY CELL	R69-0020
21	SOLENOID VALVE NO	R59-0023
22	SOLENOID VALVE NC	R59-0022
23	15 PIN CONTROLLER CORD	R69-0021
24	9 PIN CONTROLLER CORD	R69-0022

<p><b>HOSES</b> INLET P/N: 0112-0030 PRODUCT P/N: 0112-0031 REJECT P/N: 0112-0032</p> 	<p><b>GARDEN HOSE FEMALE ADAPTER KIT: NUT, WASHER, .5" BARB, P/N: 0112-0046</b></p> 	<p><b>NIPPLE</b> P/N: 041610050, 0.5 X CL, PVC P/N: 046531105, 0.5 X CL, PP P/N: 044531107, 0.5 X 1.5, SS</p> 	<p><b>ELBOW</b> P/N: 041531857 0.5 FPT, PVC P/N: 047-0002 0.5 FPT, PP</p> 
<p><b>TEE, THREAD, NPT</b> P/N: 044531606, .5, SS</p> 	<p><b>BUSHING</b> P/N: 044530212 .5M X .25F, SS P/N: 046531848 .5 M X .38 F, PP P/N: 046530201 .5 M X .25 F, PP</p> 	<p><b>NIPPLE</b> P/N: 041-0003, 0.5 X 4, PVC</p> 	<p><b>PLUG</b> P/N:10520898 3/8 STEM P/N:10521106 1/2 STEM</p> 
<p><b>BARB CONNECTORS</b> .38 STEM x .38 BARB P/N: 10521026 .5 STEM X .5 BARB P/N: 10521011</p> 	<p><b>CHECK VALVE,.5MPT</b> P/N:10720102</p> 	<p><b>CON, MALE, T X NPT, SS</b> P/N: 10-0020 .25 T X .25 NPT P/N: 10-0021 .38 T X .5 NPT</p> 	<p><b>ELBOW, MALE, T X NPT, SS</b> P/N: 10-0023 .5 T X .5 NPT</p> 
<p><b>REPLACEMENT NUTS, SS</b> 10-0027 .38" NUT 10-0028 .5" NUT</p> 	<p><b>REPLACEMENT FERRULES, SS</b> 10-0030 .38" FERRULE SET 10-0031 .5" FERRULE SET</p> <p>FITTING SIDE  NUT SIDE </p> <p>SOLD IN SETS OF 10 ONLY</p>	<p><b>MALE QUICK CONNECTORS</b> P/N: 10-L004 .38 T X .38 NPT P/N: 10-L006 .5 T X .5 NPT</p> 	<p><b>FEMALE QUICK CONNECTOR</b> P/N: 10-L026 .25 T X .25 NPT P/N: 10-L027 .38 T X .38 NPT</p> 
<p><b>SWIVEL EL, TUBE X NPT</b> P/N: 10-L107 .5 T X .5 NPT</p> 	<p><b>SWIVEL BR TEE, TUBE X NPT</b> P/N: 10-L127 .25 T X .25 NPT P/N: 10-L129 .38 T X .25 NPT</p> 	<p><b>SWIVEL RUN TEE, TUBE X NPT</b> P/N: 10-L152 .25 T X .25 NPT P/N: 10-L157 .5 T X .5 NPT</p> 	<p><b>UNION ELBOW</b> P/N: 10-L202 .38 T X .38 T</p> 
<p><b>UNION</b> P/N: 10-L230 .5 T X .5 T</p> 	<p><b>UNION TEE</b> P/N: 10-L253 .38 T X .25 T P/N: 10-L254 .5 T X .5 T P/N: 10-L255 .5 T X .38 T</p> 	<p><b>REDUCER</b> P/N: 10-L352 3/8" X 1/4" P/N: 10-L353 1/2" X 3/8"</p> 	<p><b>PLUG-IN EL, STEM x TUBE</b> P/N: 10-L376 .25 P/N: 10-L377 .38 P/N: 10-L378 .5</p> 
<p><b>PLUG-IN TEE, STEM x TUBE</b> P/N: 10-L402 .38 T X .38 STEM</p> 	<p><b>Q-CON</b> P/N: 16-0067 .38 F X .38 NPT</p> 	<p><b>FILTER WRENCH</b> P/N: 21675184</p> 	<p><b>EJECTOR</b> P/N: 40540-4</p> 
<p><b>CHECK VALVE, TUBE X TUBE</b> P/N:55-0011 .38 P/N: 55-0015 .25</p> 	<p><b>STRAIN RELIEF</b> P/N:66510619</p> 	<p><b>STRAIN RELIEF FITTING P/N:66510621</b> 3 HOLE BUSHING PN- 66-0041</p> 	<p><b>CASTER 100MM</b> P/N: 94-0013 W/BRAKE P/N: 94-0014 W/O BRAKE</p> 

## **ROUTINE REPLACEMENT ITEMS (NON-DURABLE COMPONENTS)**

<b>PART #</b>	<b>DESCRIPTION</b>
20-5102	10 Micron Carbon Block Filter
20-0009	Nephros DSU 510K Filter
24-0026	Membrane Cap, Clamp and O-ring Kit
95810125	PT401 Antiscalant, (4) 1 Gallon Containers
95810126	PT401 Antiscalant for High pH WATER, (4) 1 Gallon Containers
R22-4026	MRO1 Membrane
00810118R	MRO1 Catalytic Carbon Exchange Tank
95-0006	Peracidin (2) Quart Containers
95-0007	Peracidin (4) Quart Containers
000-095-0001	Kit,Biotrol+,5 Gal Container+Powder,Portable MRO
97HP20401	Test Strips Peracid Test (6 Bottles of 100 Strips each) For Measuring High Range Paracetic Acid
97PX20501	Test Strips Renal Check (6 Bottles of 100 Strips each) For Measuring Residual Peroxide
97PH20901	Test Strips pH (6 Bottles of 100 Strips each) For Measuring pH/Water/Acid,Base/Bicarbonate/Dialysate
97RC22101	WaterCheck RC (6 Bottles of 100 Strips each) For Measuring Residual Chlorine
97CM20201	WaterCheck 2 (6 Bottles of 100 Strips each) For Measuring Low Level Chlorine/Chloramine
LAL	LAL Endotoxin Testing, Exact results in just a few days
AAMI	AAMI Chemical Analysis, Results within one week

\*Call AmeriWater or your AmeriWater distributor for pricing.

## **SECTION 11, DISPOSAL**

Disposal of the device is the responsibility of the Medical Director of the facility. All local codes and regulations regarding the disposal must be followed.

Ameriwater recommends that the entire device be cleaned / decontaminated prior to beginning the disposal process. Many of the major components may be eligible for recycling in your area, except for the controller and membranes. It is recommended that these items be incinerated.