

ENCLOSURE BILL OF MATERIALS

1-22	66-0155	1	CONNECTOR, RECEPTACLE, 8-PIN, FEMALE, M12, PANEL MOUNT
1-21	66510621	1	CONNECTOR, CORD GRIP, 0.5NPT
1-20	64-0034	1	JUMPER, RELAY, GREY
1-19	66932111	1	LINE,CORD,120V,14/3,HOSP GRADE,STRIPPED ENDS
1-18	66-0135	4	CONNECTOR, RECEPTACLE, 4-PIN, FEMALE, M12, PANEL MOUNT
1-17	64760151	8"	DIN RAIL, 35MM
1-16	64-0033	5	RELAY, 1PDT, 24VDC, SPRING CAGE
1-15	66-0133	1	TERMINAL MARKERS, BLANK
1-14	66-0161	1	JUMPER, TERMINAL, 2 POSITION
1-13	66-0131	1	END PLATE, GREY
1-12	66-0130	2	TERMINAL BLOCK, GREEN/YELLOW, 2 TIER, 4 COND.
1-11	66-0129	2	TERMINAL BLOCK, GREY, 2 TIER, 2 COND.
1-10	66-0128	3	TERMINAL BLOCK, GREY, 2 TIER, 4 COND.
1-9	66-0127	1	END PLATE, ORANGE
1-8	63-0023	2	FUSE HOLDER, GREY
1-7	63-0022	2	END BARRIER, GREY
1-6	63-0026	1	FUSE, 2.5A, SLOW BLOW
1-5	63-0021	1	FUSE, 3A
1-4	77-0021	1	POWER SUPPLY, 24VDC, 120W
1-3	60-0053	1	BRACKET, WALL MOUNT
1-2	60-0052	1	SUBPANEL, PLATED STEEL, 8.78x5.83
1-1	60-0051	1	ENCLOSURE, POLYCARBONATE, OPAQUE, 9.61x6.65x5.59
ITEM #	AMERI WATER P/N	1	DESCRIPTION

5. REPLACEMENT PARTS LIST

Part #	Description
53-0035	HEATSAN VALVE ISOLATION INTERFACE
066-0001	CABLE, INTERFACE ISO
66-0152	CABLE, VALVE, DIN TO MENCOM
000-HVIS	HEATSAN VALVE ISOLATION SYSTEM
001-144-0001	VAL ASSY, HEATSAN ISOLATION MANIFOLD
001-144-0002	MANIFOLD ASSY, ACID SAMPLE, ISOLATION
001-144-0003	MANIFOLD ASSY, ACID SAMPLE PORT
001-144-0004	VAL ASSY, ACID ISOLATION
001-144-0005	VAL ASSY, BICARB ISOLATION
59-0054	VAL,SOL,.75NPT,NO,24VDC COIL