

## 9 UNIQUE COMPONENT PART NUMBERS

### UNIQUE COMPONENT PARTS (NON-DURABLE COMPONENTS)

PART #	DESCRIPTION
20-6001	FILTER CARTRIDGE,MIXED BED DI,20"x2.5",11000 mg/L
20-3046	FILTER CARTRIDGE ABSOLUTE ENDOTOXIN PALL,0.2 MICRON, 2.5" x 20", 222
20-5101	FILTER CARTRIDGE CARBON BLOCK, 10 MICRON, 2.5" X 10", DOE
0185-0070	TANK,VENT,FILTER,ASSY,.5",2.5x10,.75
20-3021	VENT,FILTER CARTRIDGE SUBMICRON, 0.2 MICRON, 2.5" X 10", DOE
R22-2521	MEMBRANE, 2.5 X 21 ( System requires quantity of 3 )
24-0031	O-RING,LARGE FOR 24-0024 2.5" SST HOUSING ( RO Membrane )
24-0035	O-RING, 4x26, 24-0024 END CAP 2.25, and 40' SS HSG, SMALL BLUE O-RING ( RO Membrane )
21-0035	O-RING,FLTR,HOUS,#10 "SLIMLINE", BLUE/BLK, .50I/O,NO PR ( Prefilter, Vent Filter and DI )
21530233	O-RING,FLTR HSG,10IN BB,SQ THD ( Loop Endotoxin Filter Assembly)

\*Call AmeriWater or your AmeriWater distributor for pricing.

For 00SRO1095 RO, refer to Manual 98-0150

For 00MW200 Silex Wall Mount Deionizer, refer to Manual 98-0160

For 0185-8000 or 00HC-0090-PRV Healthcare Storage Tanks, refer to Manual 098-0001

For 0021003 Loop Endotoxin Filter, refer to Manual 98-0033