



# Sterile Processing Water Systems

Upgrade your current system for the AAMI ST108 Standard

High quality water is imperative for patient health and equipment life.

Meeting the AAMI ST108 Standards means your facility is committed to preventing the spread of bacteria and other contaminants that may harm your patients or wear your equipment prematurely.

Upgrading your reliable AmeriWater system is simple and cost-efficient. Contact us today for a free consultation at 800 535 5585 or [Info@AmeriWater.com](mailto:Info@AmeriWater.com).

[www.ameriwater.com](http://www.ameriwater.com)

# Create critical water consistently with AmeriWater's AAMI ST108 Standard Upgrade

## Your upgrade will include:



### Ultraviolet Light

- Controls bacteria proliferation in purified water storage and distribution systems
- Equipped with an online monitor to measure output intensity and alarms when lamp needs replaced



### Ultrafiltration

- 20-3046 Endotoxin cartridge (included)
- Installed post UV to remove endotoxins from the purified water storage and distribution system



### Alarm Panel

- Monitors system operating conditions
- Includes RO alarm, low storage tank, and DI resistivity
- Remote alarm included for continuous system monitoring outside water room



### Storage Tank with Sealed Lid, Flowmeter and Submicron Vent Filter

Our system options come in 30 and 70 GPM. Contact our water professionals and see how this simple upgrade can save you money today.