



# AmeriWater

The Water Purification Specialists

## Sterile Processing

---

### ST108 Upgrade Kits for AmeriWater Filter Systems

- Upgrade kit includes:
  - Alarm panel
  - Tank flow meter
  - Tank vent filter
  - Ultraviolet light
  - Endotoxin filter



[www.amerewater.com](http://www.amerewater.com)

## Upgrade Your System for the AAMI ST108 Standard Simply and Cost Effectively

Having high quality water is imperative for patient health and equipment life. AmeriWater's Sterile Processing Field Upgrade Kits transform your existing AmeriWater system into one which conforms with the new AAMI ST108 Standard Water for the Processing of Medical Devices. Meeting the AAMI ST108 Standards means your facility is committed to preventing the spread of bacteria and other contaminants that may harm your patients or wear your equipment prematurely.

There are two kits available 001-121-0020 (30 gpm at 30 psi) and 001-121-0021 (70 gpm at 30 psi).

### Upgrade Kit Includes:

- Alarm Panel
  - Monitors system operating conditions
  - Includes RO alarm, low storage tank, and DI resistivity
  - Remote alarm included for continuous system monitoring outside water room
- Storage tank distribution loop return flowmeter assembly
- Storage tank submicron vent filter assembly
- Ultraviolet Light
  - Controls bacteria proliferation in purified water storage and distribution systems
  - Equipped with an online monitor to measure output intensity and alarms when lamp needs replaced
- Ultrafiltration (endotoxin cartridges included)
  - Installed post UV to remove endotoxins from the purified water storage and distribution system

DI water is required for the system in order to meet the AAMI standards, see choices below:

- Two column Silex Deionizer and DI bypass header
- DI Exchange Tanks handled by a local DI exchange provider and DI bypass header

